

Patients come in all sizes – from morbidly obese to small child, and any weight in between – that’s why there’s a complete line of GlideScope® video laryngoscopes.

The GlideScope® GVL® consistently provides a clear view, enabling quick intubations.

Proven, practical, and reliable, the GlideScope® GVL® achieves a Cormack-Lehane Grade I or Grade II view 99% of the time.¹

Whether it’s in the ED, the OR, or the ICU, patients present in all shapes and sizes – it’s important to choose an airway device that’s right for them. The GlideScope® GVL®.




Actual airway view

FEATURES

- High-resolution camera for real-time view of the airway and tube placement
- Anti-fogging mechanism to resist lens contamination
- Non-glare color monitor with NTSC video output to facilitate recording
- Unique blade angulation
- Available in 4 sizes (small child to bariatric patients)
- Operational in seconds
- Easy to use, easy to teach



**A consistently clear view
enabling quick intubations.**



REFERENCE: 1. Cooper, RM. Cardiothoracic Anesthesia, Respiration and Airway; Early clinical experience with a new video laryngoscope (GlideScope®) in 728 patients. *Canadian Journal of Anesthesia* 2005; 52: 2: 191-198.

	® 2	GVL® 3	GVL® 4	GVL® 5
Patient Weight	1.8 - 10 kg	10 kg - Adult	40 kg - Morbidly Obese	40 kg - Morbidly Obese
GlideRite® Rigid Stylet			For use with endotracheal tubes 6.0 mm and larger.	