

The GlideRite® Rigid Stylet is specifically designed to work with GlideScope® Video Laryngoscopes.

The angle of the GlideRite® Rigid Stylet complements the unique GlideScope® angle to help facilitate quick placement of an endotracheal tube and to help reduce patient trauma.

Convenient and cost-effective, the reusable GlideRite® Rigid Stylet is recommended as a valuable element of any intubation with a GlideScope® Video Laryngoscope.



## FEATURES

- ▶ Provides improved maneuverability for placement of endotracheal tube
- ▶ Complements unique GlideScope® angle
- ▶ Constructed of durable, reusable stainless steel
- ▶ Easy to clean or sterilize in an autoclave
- ▶ Can be used for up to 250 autoclave cycles, or one year of use
- ▶ Designed for use in ET tubes 6mm and larger
- ▶ Available in convenient packs of 6
- ▶ Easy to learn, use & teach



## Specifically designed to work with GlideScope® Video Laryngoscopes

- Stylet/ETT insertion can be viewed real-time via GlideScope® system's integrated high-resolution camera and monitor, helping ensure accurate ETT placement
- GlideScope® Video Laryngoscopes offer a consistently clear view of the airway, enabling quick intubation
- Reusable and single-use GlideScope® configurations available
- Available in a variety of sizes, for patients from morbidly obese to preterm, and any weight in between

GlideScope®, GlideRite®, GVL®, Verathon® and Verathon Medical® are registered trademarks of Verathon Inc. in the USA and/or other countries. All rights reserved. Copyright © 2008 by Verathon Inc.

GlideScope® video laryngoscope systems are CE marked in accordance with the Medical Device Directive, and the Verathon Inc. quality system is Quality System Certified to ISO 13485:2003 standards. US Patent No. 6,655,377; 6,543,447; International patents pending. 0900-1281-03-86