


ScanPoint® Bladder Imaging Technology

ONLINE and LOCAL CLIENT Versions



ScanPoint® Online

(For BladderScan® bladder volume (BVI) and mass (BVM) instruments - BVI 6100, BVI 6200, BVI 6300, BVI 6400 & BVM 6500)

Scan
Point®

Verathon™ designed ScanPoint®, an innovative, online service, to enhance and expand the capabilities of the BladderScan® 6000 series of portable ultrasound bladder volume and mass measurement devices.

Using ScanPoint® Online you can:

- View ultrasound measurements and images from patient exams
- Print exam results for patient records and for reimbursement (if applicable) using your choice of report formats
- Save and store images online and take advantage of the ScanPoint® automatic file backup service
- Ensure the optimum accuracy of your instrument by calibrating online, conveniently, in a matter of minutes and as often as you wish
- Enjoy software upgrades automatically
- Make use of Verathon™ programs including:
 - Warranty against scanner defects
 - Insurance against damage, theft and loss (deductible applies)
 - Replacement with a new unit in the event of accidental droppage (deductible applies)



ScanPoint® Total ReliabilitySM Plan

The ScanPoint® Plan includes:

- ScanPoint® online anytime calibration and certification service for your BladderScan® BVI 6000 series instrument
- Technology assurance
- Software assurance
- Instrument insurance
- Coverage of all shipping charges

With the ScanPoint® Plan, you can calibrate your instrument online for optimum accuracy. Simply scan the calibration targeting system and use the communication cradle to transmit calibration data over the Internet.


While under continuous coverage from the ScanPoint® program, your scanning unit is warranted perpetually against outdated technology, which means we'll send you a new unit without charge when your instrument model is no longer manufactured and cannot be repaired. You'll also receive software updates for your BladderScan®, which you can download from the ScanPoint® Website. That way you'll always have the latest version of the ScanPoint® System for Bladder Volume Measurement.

Noninvasive, Accurate, Reliable, & Easy to Use



ScanPoint® Bladder Imaging Technology

ONLINE and LOCAL CLIENT Versions



ScanPoint® Local Client (For BladderScan® bladder volume (BVI) instruments - BVI 6100, BVI 6200 & BVI 6400)



If you are not ready to take advantage of ScanPoint® Online bladder imaging, consider **ScanPoint® Local Client (LC)**. With ScanPoint® LC:

- You can view, print and save images on your desktop.
- Calibration is included.
- A warranty against scanner defects is offered, as well as insurance against damage, theft and loss, and replacement with a new unit in the event of accidental droppage.

Benefit	ScanPoint® LC (Local Client)	ScanPoint® Online
Choosing the ScanPoint® service that best suits your practice:	For the practice that: <ul style="list-style-type: none"> • Prefers access to exam records using their patient data management program, systems and reporting templates • May not have internet access to take advantage of ScanPoint® Online imaging services • Doesn't need the immediacy of online software updates or calibration. Can accommodate unit being shipped to VM for calibration service • Doesn't require imaging services for the BVI 6300 and/or the new BVM 6500 instruments 	For the practice that: <ul style="list-style-type: none"> • Prefers the convenience and speed of Internet-based imaging services • Prefers the less-hardware-intensive approach of online exam record storage • Prefers the immediacy and frequency of both calibration online, and device support • Prescribes the BVI 6300 in addition to the other BladderScan® models • Uses the new BVM 6500, to measure ultrasound-estimated bladder weight (UEBW)
Note:	Requires customer capability to ensure security and integrity of exam record storage Requires routine data backup	Automatically offers patient privacy protections Data backup is routine and automatically provided
View images	BVI 6100, 6200, 6400	BVI 6100, 6200, 6400 Plus BVI 6300 & BVM 6500
Print images	BVI 6100, 6200, 6400	BVI 6100, 6200, 6400 Plus BVI 6300 & BVM 6500
Save images	BVI 6100, 6200, 6400	BVI 6100, 6200, 6400 Plus BVI 6300 & BVM 6500
Storage of data (exam records)	Via Local Client in local format Customer must make additional storage requirements	Via Internet VM provides offsite storage at no additional charge
Calibration	Available	Instant via the Internet Calibration can be as often as desired VM has ability to resolve issues real time
Software	Install CD provided	"Thin client" – No additional software required
ScanPoint® upgrades	Available via mailed CD	Easily accessible online
Scanner upgrades	Handled when the unit is sent to VM for calibration	Updated as soon as scanner connects to ScanPoint® Online

BladderScan®, ScanPoint®, Verathon™ and Verathon Medical™ are either registered trademarks or trademarks of Verathon Inc. in the USA and/or other countries. Copyright ©2006. All rights reserved.



Verathon Medical (United Kingdom) Ltd.

The Granary Manor Farm Courtyard / Aston Sandford, Aylesbury
Buckinghamshire, HP17 8JB / United Kingdom
+44.1844.299.207 / +44.1844.299.218 Fax
www.verathon.com